
**Instructions for use**  
**2-MET Plasma RIA** **Fast Track**

Please use only the valid version of the Instructions for Use provided with the kit

**REF****BA R-8300****IVD****400 kBq**

## 1. Introduction

### 1.1 Intended use and principle of the test

<sup>125</sup>I – Radioimmunoassay for the quantitative determination of free Metanephrine and free Normetanephrine in plasma.

#### **Related Products:**

2-MET Plasma ELISA <sup>Fast Track</sup>	
Metanephrine Plasma ELISA <sup>Fast Track</sup>	Metanephrine Plasma RIA <sup>Fast Track</sup>
Normetanephrine Plasma ELISA <sup>Fast Track</sup>	Normetanephrine Plasma RIA <sup>Fast Track</sup>

Alternatively the parameters can also be run automatically on an ELISA processor such as the Gemini instrument from Stratec Biomedical. The Gemini protocol is available upon request.

Metanephrine (Metadrenaline) and Normetanephrine (Normetadrenaline) are first extracted using an ion exchange matrix followed by an acylation process.

The assay procedure follows the basic principle of radioimmunoassay, involving competition between a radioactive and a non-radioactive antigen for a fixed number of antibody binding sites. The amount of <sup>125</sup>I-labelled antigen bound to the antibody is inversely proportional to the analyte concentration of the sample. When the system is in equilibrium, the antibody bound radioactivity is precipitated with a second antibody in the presence of polyethylene glycol. The precipitate is counted in a gamma counter. Quantification of unknown samples is achieved by comparing their activity with a reference curve prepared with known standards.

⚠ *The antibodies used in this test kit only recognise the biologically relevant L-forms of Metanephrines. Commercially available synthetic Normetanephrine or Metanephrine is always a mixture of the D- and L-form. The ratio between both forms differs widely from lot to lot. This has important implications if synthetic Metanephrines are used to enrich native samples. As only about 50% of the synthetic Metanephrines - the L-portion - will be detected by use of this kit, spiked samples will be underestimated. Therefore native samples containing solely the L-form should be used.*

### 1.2 Clinical application

Metanephrine and Normetanephrine are the metabolites of the catecholamines Epinephrine and Norepinephrine, respectively. Cells derived from neuroendocrine tumors (e.g. pheochromocytoma) are known to produce catecholamines, which are secreted episodically via vesicles into the blood stream. But beside this, a small portion of the catecholamines is metabolized inside the cells to the corresponding catecholamines metabolites – namely Metanephrine, Normetanephrine and 3-Methoxytyramine – which are secreted at low levels continuously into the blood stream.

Recent studies and publications have shown that the quantification of these plasma free Metanephrine and plasma free Normetanephrine is the most accurate biochemical marker for the clinical diagnosis of pheochromocytoma and follow-up of pheochromocytoma patients.

Therapeutic consequences should never be based on laboratory results alone, even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient.

Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient, it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

## 2. Procedural cautions, guidelines, warnings and limitations

### 2.1 Precautions, guidelines and warnings

- (1) This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- (2) This assay was validated for a certain type of sample as indicated in *Intended Use* (please refer to Chapter 1). Any off-label use of this kit is in the responsibility of the user and the manufacturer cannot be held liable.
- (3) Reagents of this kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
- (4) The principles of Good Laboratory Practice (GLP) have to be followed.
- (5) In order to reduce exposure to potentially harmful substances, wear lab coats, disposable protective gloves and protective glasses where necessary.

- (6) All kit reagents - with the exception of Precipitating Reagent – and specimens should be brought to room temperature and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- (7) For dilution or reconstitution purposes use deionized, distilled, or ultra-pure water.
- (8) The radioactive material (<sup>125</sup>Iodine, half life 60 days, emitting ionizing X-radiation with 28 keV and G-radiation with 35.5 keV) may be received, acquired, possessed and used only by physicians, laboratories or hospitals. Products are dispatched on the basis of the nuclear and radiation protection regulations.
- (9) Duplicate determination of sample is highly recommended to be able to identify potential pipetting errors.
- (10) Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- (11) Incubation times do influence the results. All tubes should be handled in the same order and time intervals.
- (12) To avoid cross-contamination of reagents, use new disposable pipette tips for dispensing each reagent, sample, standard and control.
- (13) A standard curve must be established for each run.
- (14) The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report.
- (15) Do not mix kit components with different lot numbers within a test and do not use reagents beyond expiry date as shown on the kit labels.
- (16) For information on hazardous substances included in the kit please refer to Safety Data Sheet (SDS). The Safety Data Sheet for this product is available directly on the website of the manufacturer or upon request.
- (17) The expected reference values reported in this test instruction are only indicative. It is recommended that each laboratory establishes its own reference intervals.
- (18) The results obtained with this test kit should not be taken as the sole reason for any therapeutic consequence (e.g. medication before a scheduled surgery) but have to be correlated to other diagnostic tests and clinical observations.
- (19) Kit reagents must be regarded as hazardous waste and disposed according to national regulations.

## 2.2 Limitations

Any inappropriate handling of samples or modification of this test might influence the results.

### 2.2.1 Interfering substances

Samples containing precipitates or fibrin strands or which are haemolytic or lipemic might cause inaccurate results.

### 2.2.2 Drug interferences

Please refer to point "Sample collection and storage".

### 2.2.3 High-Dose-Hook effect

No hook effect was observed in this test.

## 3. Storage and stability

Store the unopened reagents at 2 - 8 °C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened, the reagents are stable for 1 month when stored at 2 - 8 °C.

## 4. Materials

### 4.1 Content of the kit

<b>BA D-0090</b>	<b>FOILS</b>	<b>Adhesive Foil</b> - Ready to use
Contents:	Adhesive Foils in a resealable pouch	
Volume:	1 x 4 foils	
<b>BA R-0030</b>	<b>PREC-REAG</b>	<b>Precipitating Reagent</b> - Ready to use
Content:	Goat anti-rabbit serum in PEG phosphate buffer	
Volume:	2 x 55 ml/vial, yellow cap	

**BA R-0120** <sup>125</sup>I **ADR MN** <sup>125</sup>I – **Metanephrine** - Ready to use  
 Content: <sup>125</sup>I labeled Metanephrine, red coloured  
 Volume: 1 x 5.5 ml/vial, blue cap  
 Hazards identification:  Radioactive, activity < 200 kBq

**BA R-0220** <sup>125</sup>I **NAD NMN** <sup>125</sup>I – **Normetanephrine** - Ready to use  
 Content: <sup>125</sup>I labeled Normetanephrine, red coloured  
 Volume: 1 x 5.5 ml/vial, yellow cap  
 Hazards identification:  Radioactive, activity < 200 kBq

**Standards and Controls** - Ready to use

Cat. no.	Component	Colour/Cap	Concentration pg/ml		Concentration pmol/l		Volume/Vial
			MN	NMN	MN	NMN	
<b>BA E-8301</b>	STANDARD A	white	0	0	0	0	4 ml
<b>BA E-8302</b>	STANDARD B	light yellow	36	72	183	393	4 ml
<b>BA E-8303</b>	STANDARD C	orange	120	240	608	1 310	4 ml
<b>BA E-8304</b>	STANDARD D	dark blue	360	720	1 825	3 931	4 ml
<b>BA E-8305</b>	STANDARD E	light grey	1 200	2 400	6 084	13 104	4 ml
<b>BA E-8306</b>	STANDARD F	black	3 600	7 200	18 252	39 312	4 ml
<b>BA E-8351</b>	CONTROL 1	light green	Refer to QC report for expected value and acceptable range!				4 ml
<b>BA E-8352</b>	CONTROL 2	dark red	Refer to QC report for expected value and acceptable range!				4 ml

Conversion: Metanephrine (pg/ml) x 5.07 = Metanephrine (pmol/l)  
 Normetanephrine (pg/ml) x 5.46 = Normetanephrine (pmol/l)  
 Content: Acidic buffer with non-mercury stabilizer, spiked with a defined quantity of Metanephrine and Normetanephrine

**BA R-8110** **AS MN** **Metanephrine Antiserum** - Ready to use  
 Content: Rabbit anti- metanephrine antibody, blue coloured  
 Volume: 1 x 5.25 ml/vial, blue cap

**BA R-8210** **AS NMN** **Normetanephrine Antiserum** - Ready to use  
 Content: Rabbit anti- normetanephrine antibody, yellow coloured  
 Volume: 1 x 5.25 ml/vial, yellow cap

**BA R-8312** **ACYL-CONC** **Acylation Concentrate** – Concentrated  
 Content: Acylation reagent in DMSO  
 Volume: 1 x 1.5 ml/vial, dark grey cap

Hazards identification:  H302 Harmful if swallowed.

**BA R-8313** **ASSAY-BUFF** **Assay Buffer** - Ready to use  
 Content: 25 % organic solvent  
 Volume: 1 x 30 ml/vial, orange cap

**BA R-8318** **EXTRACT-PLATE 96** **Extraction Plate** - Ready to use  
 Content: 1 x 96 well (6x16) plate, precoated with ion-exchanger in a resealable pouch

**BA R-8325** **CLEAN-CONC 25x** **Cleaning Concentrate** - Concentrated 25x  
 Content: Buffer with sodium acetate  
 Volume: 1 x 20 ml/vial, brown cap

**BA R-8326**    **ELUTION-BUFF**    **Elution Buffer** - Ready to use

Content: 0.1 M Sodium hydroxide, dark purple coloured

Volume: 1 x 14 ml/vial, dark green cap

**BA R-8327**    **ADJUST-BUFF**    **Adjustment Buffer** - Ready to use

Content: TRIS buffer with BSA and a non-mercury stabilizer

Volume: 1 x 6 ml/vial, light purple cap

**BA R-8828**    **EQUA-REAG**    **Equalizing-Reagent** - Ready to use

Content: Human serum, negative for HIV I/II, HBsAg and HCV

Volume: 1 x 14 ml/vial, white cap

#### 4.2 Additional materials and equipment required but not provided in the kit

- Calibrated precision pipettes to dispense volumes between 20 - 500 µl; 3 ml
- Conical tubes and suitable rack
- Centrifuge capable of at least 3 000 x g
- Suitable device for aspirating or decanting the tubes
- Microtiter plate shaker (shaking amplitude 3 mm; approx. 600 rpm)
- Gamma Counter
- Vortex mixer
- Water (deionized, distilled, or ultra-pure)

For the washing steps during the extraction a microtiter plate washing device (manual, semi-automated or automated) may be applied.

#### 5. Sample collection and storage

Medications like Serotonin-noradrenaline reuptake inhibitors, tricyclic antidepressants, MAO inhibitors, antihypertensive drugs and L-DOPA can influence Metanephrine and Normetanephrine level. People who are taking such medication should consult with their doctor before specimen collection.

Sympathomimetic agents, sport and smoking can also influence Metanephrine and Normetanephrine level.

Alcohol and caffeinated drinks should be avoided the day before and including the day of sample collection.

##### **EDTA- or Heparin-Plasma**

Whole blood should be collected into centrifuge tubes (Monovette™ or Vacuette™) containing EDTA or heparin as anti-coagulant and centrifuged (according to manufacturer's instructions) immediately after collection.

Haemolytic and lipemic samples should not be used for the assay.

Storage: up to 6 hours at 2 - 8 °C, for longer period (up to 6 months) at -20 °C.

Repeated freezing and thawing should be avoided.

#### 6. Test procedure

Allow all reagents - with the exception of Precipitating Reagent - to reach room temperature and mix thoroughly by gentle inversion before use. Number the extraction wells / RIA tubes accordingly. Duplicates are recommended.

 *Pipetted liquids should not adhere to the wall of the RIA tubes. If necessary please centrifuge the tubes for 1 minute at 500 x g to spin down adhering liquids.*

##### 6.1 Preparation of reagents

###### **Cleaning Buffer**

Dilute the 20 ml Cleaning Concentrate with water (deionized, distilled, or ultra-pure) to a final volume of 500 ml.

Storage: 1 month 2 - 8 °C

###### **Acylation Solution**

As the Acylation Solution is only **stable for a maximum of 3 minutes**, it should not be prepared before starting the assay. Therefore its preparation is described in the protocol in chapter 6.3, step 4 and chapter 6.4, step 4.

Discard after use!

## 6.2 Preparation of samples

 The extraction procedure is the same for Metanephrine and Normetanephrine and has to be done only once.

### Extraction

The following extraction procedure can be run with 200 µl or 250 µl of plasma sample.

The procedure for 250 µl plasma is highlighted in grey and italicised and may be used in case higher supernatant volumes for pipetting to the subsequent RIA are preferred.

The RIA procedure itself is not affected by this alternative protocol.

1.	Pipette <b>20 µl</b> of <b>standards</b> and <b>controls</b> into the respective wells of the <b>Extraction Plate</b> . <b>Alternatively</b> pipette <b>25 µl</b> of <b>standards</b> and <b>controls</b> into the respective wells of the <b>Extraction Plate</b> .		
2.	Add <b>20 µl Standard A</b> to all wells containing <b>plasma samples</b> . <b>Alternatively</b> add <b>25 µl Standard A</b> to all wells containing <b>plasma samples</b> .		
3.	Add <b>200 µl</b> of <b>Equalizing Reagent</b> to the wells with <b>standards and controls</b> . <b>Alternatively</b> add <b>250 µl</b> of <b>Equalizing Reagent</b> to the wells with <b>standards and controls</b> .		
4.	Pipette <b>200 µl</b> of <b>plasma samples</b> to the respective wells. <b>Alternatively</b> pipette <b>250 µl</b> of <b>plasma samples</b> to the respective wells.		
5.	Incubate plate for <b>2 hours</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm)		
6.	Empty plate and blot dry by tapping the inverted plate on absorbent material		
7.	Pipette <b>250 µl</b> of <b>Assay Buffer</b> into all wells. Incubate the plate for <b>5 min</b> at <b>RT</b> (20 – 25 °C) on a <b>shaker</b> (approx. 600 rpm). Empty plate and blot dry by tapping the inverted plate on absorbent material.		
8.	Wash the plate <b>3 x</b> by adding <b>350 µl</b> of <b>Cleaning Buffer</b> , <b>discarding</b> the content and <b>blotting dry each time</b> by tapping the inverted plate on absorbent material.		
9.	Pipette <b>100 µl</b> of <b>Elution Buffer</b> into all wells. <b>Alternatively</b> pipette <b>125 µl</b> of <b>Elution Buffer</b> into all wells. <i>Please note: the colour changes caused by the elution buffer can vary between standards and samples.</i>		
10.	Cover plate with adhesive foil. Incubate <b>15 min</b> at <b>RT</b> (20 - 25 °C) on a <b>shaker</b> (approx. 600 rpm). Remove the foil. <b>Do not decant the supernatant thereafter!</b>		
	The following volumes of the supernatant are needed for the subsequent RIA: <table border="1" style="display: inline-table; vertical-align: middle;"> <tr> <td><b>Metanephrine 60 µl</b></td> <td><b>Normetanephrine 20 µl</b></td> </tr> </table>	<b>Metanephrine 60 µl</b>	<b>Normetanephrine 20 µl</b>
<b>Metanephrine 60 µl</b>	<b>Normetanephrine 20 µl</b>		

## 6.3 Metanephrine RIA

 The use of conical tubes for the RIA is highly recommended!

1.	Pipette <b>60 µl</b> of <b>Elution Buffer</b> into the <b>tubes</b> for the <b>NSB</b> .
2.	Pipette <b>60 µl</b> of the <b>extracted standards, controls</b> and <b>samples</b> into the respective tubes. <i>Please hold the Extraction Plate at a slight angle in order to facilitate this pipetting step.</i>
3.	Pipette <b>25 µl</b> of <b>Adjustment Buffer</b> into <b>all tubes (except totals)</b> .
4.	Preparation of <b>Acylation Solution</b> : Pipette <b>80 µl Acylation Reagent Concentrate</b> (BA R-8312) to <b>3 ml water</b> (deionized, distilled, or ultra-pure) and mix thoroughly.
5.	Pipette <b>25 µl</b> of the freshly prepared <b>Acylation Solution</b> into <b>all tubes (except totals)</b> .
6.	Mix thoroughly (vortex) and incubate for <b>15 min</b> at <b>RT</b> (20 - 25 °C).
7.	Pipette <b>50 µl</b> of <b>Metanephrine Antiserum</b> into <b>all tubes (except totals and NSB)</b> ; mix thoroughly (vortex).
8.	Incubate for <b>1 h</b> at <b>RT</b> (20 - 25 °C).
9.	Pipette <b>50 µl</b> of the <b><sup>125</sup>I Metanephrine</b> into <b>all tubes</b> and mix thoroughly (vortex).
10.	Cover tubes. Incubate for <b>15 - 20 h</b> (overnight) at <b>2 - 8 °C</b> .
11.	Mix the chilled (2 - 8 °C) <b>Precipitating Reagent</b> thoroughly, pipette each <b>500 µl</b> into <b>all tubes (except totals)</b> , and mix on a vortex.
12.	Incubate for <b>15 min</b> at <b>2 - 8 °C</b> .
13.	Centrifuge for <b>15 min</b> at <b>3 000 x g</b> , if possible in a refrigerated centrifuge.
	<b>Continue without any delay with step 14.</b>
14.	<b>Decant</b> or aspirate the <b>supernatant carefully (except totals)</b> . Blot the tubes dry and leave them upside down for 2 minutes.
15.	<b>Count</b> all tubes for <b>1 min</b> in a gamma counter.

## 6.4 Normetanephrine RIA

 The use of conical tubes for the RIA is highly recommended!

1.	Pipette <b>20 µl</b> of <b>Elution Buffer</b> into the <b>tubes</b> for the <b>NSB</b> .
2.	Pipette <b>20 µl</b> of the <b>extracted standards, controls</b> and <b>samples</b> into the respective tubes. <i>Please hold the Extraction Plate at a slight angle in order to facilitate this pipetting step.</i>
3.	Pipette <b>25 µl</b> of <b>Adjustment Buffer</b> into <b>all tubes (except totals)</b> .
4.	Preparation of <b>Acylation Solution</b> : Pipette <b>80 µl Acylation Reagent Concentrate</b> (BA R-8312) to <b>3 ml water</b> (deionized, distilled, or ultra-pure) and mix thoroughly.
5.	Pipette <b>25 µl</b> of the freshly prepared <b>Acylation Solution</b> into <b>all tubes (except totals)</b> .
6.	Mix thoroughly (vortex) and incubate for <b>15 min</b> at <b>RT</b> (20 - 25 °C).
7.	Pipette <b>50 µl</b> of <b>Normetanephrine Antiserum</b> into <b>all tubes (except totals and NSB)</b> ; mix thoroughly (vortex).
8.	Incubate for <b>1 h</b> at <b>RT</b> (20 - 25 °C).
9.	Pipette <b>50 µl</b> of the <b><sup>125</sup>I Normetanephrine</b> into <b>all tubes</b> and mix thoroughly (vortex).
10.	Cover tubes. Incubate for <b>15 - 20 h</b> (overnight) at <b>2 - 8 °C</b> .
11.	Mix the chilled (2 - 8 °C) <b>Precipitating Reagent</b> thoroughly, pipette each <b>500 µl</b> into <b>all tubes (except totals)</b> , and mix on a vortex.
12.	Incubate for <b>15 min</b> at <b>2 - 8 °C</b> .
13.	Centrifuge for <b>15 min</b> at <b>3 000 x g</b> , if possible in a refrigerated centrifuge.  <b>Continue without any delay with step 14.</b>
14.	<b>Decant</b> or aspirate the <b>supernatant carefully (except totals)</b> . Blot the tubes dry and leave them upside down for 2 minutes.
15.	<b>Count</b> all tubes for <b>1 min</b> in a gamma counter.

## 7. Calculation of results

Measuring range	Metanephrine	Normetanephrine
	6.4 - 3 600 pg/ml	24.1 - 7 200 pg/ml

Subtract the mean cpm of the non-specific binding NSB from the mean cpm of standards, controls and samples.

The standard curve, from which the concentrations in the samples can be taken, is obtained by using the percentage of (B-NSB)/(B<sub>0</sub>-NSB) measured for the standards (linear, y-axis) against the corresponding concentrations (logarithmic, x-axis).

Use non-linear regression for curve fitting (e.g. spline, 4- parameter, akima).

 This assay is a competitive assay. This means: the counts are decreasing with increasing concentrations of the analyte. Counts found below the standard curve correspond to high concentrations of the analyte in the sample and have to be reported as being positive.

The concentrations of the **samples** and **controls** can be read directly from the standard curve.

Samples found with concentrations higher than the highest standard (Standard F) should be diluted accordingly with the included Equalizing Reagent and have to be re-assayed.

### Conversion

Metanephrine (pg/ml) x 5.07 = Metanephrine (pmol/l)

Normetanephrine (pg/ml) x 5.46 = Normetanephrine (pmol/l)

### Expected reference values

It is strongly recommended that each laboratory should determine its own reference values.

The expected reference values indicated below are based on method comparison studies to LC-MS/MS <sup>(1)</sup> with blood samples taken in the sitting position.

Metanephrine	Normetanephrine
< 65 pg/ml	< 196 pg/ml

For the interpretation of the results, a grey area has to be considered. This grey area does not depend on the methodology used and is reflected in a slight to moderate increase in Metanephrine and Normetanephrine up to 4 times the upper cut-off (Eisenhofer et al. 2003). Approx. 20 % of the tumors are found in this grey area, especially in the case of the Hereditary Syndrome, incidental tumors and in sporadic cases of Pheochromocytomas with a diameter less than 1 cm.


In case of a result in the grey area, it is recommended to collect a new sample together with an anamnesis concerning especially influences like the medication and age of the patient.

## 7.1 Quality control

It is recommended to use control samples according to national regulations. Use controls at both normal and pathological levels. The kit controls, or other commercially available controls should fall within established confidence limits. The confidence limits of the kit controls are printed on the QC-Report.

## 7.2 Typical standard curves

 Examples, do not use for calculation!

## 8. Assay characteristics

Analytical Sensitivity	Metanephrine		Normetanephrine
	LOD (pg/ml)	5.8	21.4
	LOQ (pg/ml)	6.4	24.1

Analytical Specificity (Cross Reactivity)	Substance	Cross Reactivity (%)	
		Metanephrine	Normetanephrine
	Derivatized Metanephrine	100	1.3
Derivatized Normetanephrine	0.31	100	
3-Methoxytyramine	0.04	7.7*	
Adrenaline	< 0.01	< 0.01	
Noradrenaline	< 0.01	< 0.01	
Dopamine	< 0.01	< 0.01	
Vanillic mandelic acid	< 0.01	< 0.01	
Homovanillic acid	< 0.01	< 0.01	
L-DOPA	< 0.01	< 0.01	
L-Tyrosin	< 0.01	< 0.01	
Tyramine	< 0.01	< 0.01	

\*Normetanephrine concentrations are not influenced by 3-Methoxytyramine in case of normal 3-methoxytyramine concentrations. Only very high 3-Methoxytyramine concentrations found in rare cases of exclusively dopamine producing tumours can cause false positive results.

Precision							
Intra-Assay				Inter-Assay			
	Sample	Mean (pg/ml)	CV (%)		Sample	Mean (pg/ml)	CV (%)
Metanephrine	1	59.6	13.6	Metanephrine	1	57.8	15.1
	2	120	8.0		2	110	10.5
	3	354	6.2		3	334	13.1
	4	906	10.5		4	903	10.4
Normetanephrine	1	148	8.4	Normetanephrine	1	144	12.0
	2	293	5.3		2	271	9.5
	3	870	5.3		3	839	8.4
	4	2 197	7.6		4	2 124	8.3

Linearity		Serial dilution up to	Mean (%)	Range (%)
	Metanephrine	1:64	111	104 - 115
	Normetanephrine	1:64	109	106 - 119

Recovery		Mean (%)	Range (%)
	Metanephrine	97	89 - 103
	Normetanephrine	110	97 - 119


Method Comparison: RIA vs. LC-MS/MS <sup>(1)</sup>	Metanephrine	$y = 1.1x - 6.3$	$r^2 = 0.97; n = 47$
	Normetanephrine	$y = 0.93x - 8.4$	$r^2 = 0.98; n = 49$

## 9. References/Literature

- (1) De Jong et al. Plasma free metanephrine measurement using automated online solid phase extraction HPLC-Tandem mass spectrometry. Clin Chem, 53(9): 1684-1693 (2007)
- (2) Eisenhofer et al. Laboratory evaluation of pheochromocytoma and paraganglioma. Clin Chem, 60:1486-1499 (2014)
- (3) Eisenhofer et al. Plasma metadrenalines: Do they provide useful information about sympatho-adrenal function and catecholamine metabolism? Clin Sci (Lond), 88:533-542 (1995)
- (4) Berkel et al. Diagnosis of endocrine disease: Biochemical diagnosis of pheochromocytoma and paraganglioma. Eur J Endocrinol, 170: R109-R119
- (5) Manz et al. Development of enantioselective immunoassays for free plasma metanephrines. Ann.N.Y.Acad.Sci., 1018:582-587 (2004)
- (6) De Jong et al. Dietary Influences on Plasma and Urinary Metanephrines: Implications for Diagnosis of Catecholamine-Producing Tumors. J Clin Endocrinol Metab, 94(8):2841-2849 (2009)
- (7) Deutschbein et al. Influences of Various Confounding Variable and Storage Conditions on Metanephrine and Normetanephrine Levels in Plasma. Clin Endocrinol, 72(2):153-160 (2010)
- (8) Eisenhofer et al. Biochemical Diagnosis of Pheochromocytoma: How to Distinguish True- from False-Positive Test Results. The Journal of Clinical Endocrinology & Metabolism, 88(6):2656-2666 (2003)

⚠ For updated literature or any other information please contact your local supplier.

### Symbols:

	Storage temperature		Manufacturer		Contains sufficient for <n> tests
	Expiry date		Batch code		For in-vitro diagnostic use only!
	Consult instructions for use		Content		CE labelled
	Caution		Catalogue number		

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#### Artverwandte Produkte:

2-MET Plasma ELISA <sup>Fast Track</sup>	
Metanephrin Plasma ELISA <sup>Fast Track</sup>	Metanephrin Plasma RIA <sup>Fast Track</sup>
Normetanephrin Plasma ELISA <sup>Fast Track</sup>	Normetanephrin Plasma RIA <sup>Fast Track</sup>

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
## 9. Referenzen/Literatur

- (1) De Jong et al. Plasma free metanephrine measurement using automated online solid phase extraction HPLC-Tandem mass spectrometry. Clin Chem, 53(9):1684-1693 (2007)
- (2) Eisenhofer et al. Laboratory evaluation of pheochromocytoma and paraganglioma. Clin Chem, 60:1486-1499 (2014)
- (3) Eisenhofer et al. Plasma metadrenalines: Do they provide useful information about sympatho-adrenal function and catecholamine metabolism? Clin Sci (Lond), 88:533-542 (1995)
- (4) Berkel et al. Diagnosis of endocrine disease: Biochemical diagnosis of pheochromocytoma and paraganglioma. Eur J Endocrinol, 170: R109-R119
- (5) Manz et al. Development of enantioselective immunoassays for free plasma metanephrines. Ann.N.Y.Acad.Sci., 1018:582-587 (2004)
- (6) De Jong et al. Dietary Influences on Plasma and Urinary Metanephrines: Implications for Diagnosis of Catecholamine-Producing Tumors. J Clin Endocrinol Metab, 94(8):2841-2849 (2009)
- (7) Deutschbein et al. Influences of Various Confounding Variable and Storage Conditions on Metanephrine and Normetanephrine Levels in Plasma. Clin Endocrinol, 72(2):153-160 (2010)
- (8) Eisenhofer et al. Biochemical Diagnosis of Pheochromocytoma: How to Distinguish True- from False-Positive Test Results. The Journal of Clinical Endocrinology & Metabolism, 88(6):2656-2666 (2003)

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