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
12. Referenzen /Literatur

1. Tsang BK et al. Steroid biosynthesis by isolated human ovarian follicular cells in vitro. *J Clin Endocrinol Metab* 1980, 51, 1407 – 11.
2. Gore-Langton RE et al. Follicular steroidogenesis and its control. In: *The physiology of Reproduction*. 1988, Ed.: Knobil et al., 331-385, Raven press, New York.
3. Hall, PF. Steroid synthesis: Organization and Regulation. In: *The Physiology of Reproduction*. 1988, Ed.: Knobil et al., 975 – 988, Raven press, New York
4. Salazar-Hna CJ, Remage-Healey L, and Schlinger BA. "Synaptocrine signaling: steroid synthesis and action at the synapse." *Endocrine reviews* 32.4, 2011, 532–549.
5. Wu CH et al. Free and protein-bound plasma estradiol-17 beta during the menstrual cycle". *The Journal of Clinical Endocrinology and Metabolism*. 1976, 43 (2), 436–45.
6. Fox EM et al. Novel actions of estrogen to promote proliferation: integration of cytoplasmic and nuclear pathways. *Steroids*, 2009, 74(7), 622-627.
7. Fernando RI et al. Estradiol abrogates apoptosis in breast cancer cells through inactivation of BAD: Ras-dependent nongenomic pathways requiring signaling through ERK and Akt. *Mol Biol Cell*. 2004, 15(7), 3266-3284.
8. Bovenkerk S et al. Induction of p53 expression and function by estrogen in osteoblasts. *Calcif Tissue Int*. 2003, 73(3), 274-280.

9. Vazquez Rodriguez G et al. Estradiol promotes breast cancer cell migration via recruitment and activation of neutrophils. *Cancer Immunol Res.* 2017, 5(3), 234-247.
10. Linton L et al. Associations of serum levels of sex hormones in follicular and luteal phases of the menstrual cycle with breast tissue characteristics in young women. *PLoS One.* 2016, 11(10), e0163865.
11. Kim JY et al. Prognostic effect of preoperative serum estradiol level in postmenopausal breast cancer. *BMC Cancer.* 2013, 13:503.
12. Farahani H et al. Serum and Saliva Levels of Cancer Antigen 15-3, Carcinoembryonic Antigen, Estradiol, Vaspin, and Obestatin as Biomarkers for the Diagnosis of Breast Cancer in Postmenopausal Women. *Laboratory Medicine.* 2020, 51, 620-627.
13. Siiteri PK et al. The serum transport of steroid hormones. *Rec Progr Horm Res* 1982, 38, 457 – 510
14. Martin B. et al. Binding of steroids by proteins in follicular fluid of the human ovary. *J Clin Endocrinol Metab* 1982, 35, 443 – 447.
15. Celec P et al. Salivary Sex Hormones during the Menstrual Cycle. *Endocrine Journal* 2009, 56 (3), 521-523.
16. Gandara BK, Leresche L, and Mancl L. Patterns of Salivary Estradiol and Progesterone across the Menstrual Cycle. *Ann N Y Acad Sci.* 2007, 1098, 446–450.
17. Baird DT Ovarian steroid secretion and metabolism in women. In: *The Endocrine Function of the human ovary.* Eds.: James V.H.T., Serio M. and Guisti G.,1976, 125 – 33, Academic press, New York
18. McNtty KP et al. Concentration of estrogens and androgens in human ovarian venous plasma and follicular fluid throughout the menstrual cycle. *J. Endocrinol.* 1976, 71, 77 –85.
19. Gruhn JG, Kazer RR. *Hormonal Regulation of the Menstrual Cycle: The Evolution of Concepts.* Springer Science & Business Media.2013, 69–73. ISBN 978-1-4899-3496-3.
20. March CM et al. Roles of estradiol and progesterone in eliciting mid-cycle luteinizing hormone and follicle stimulating hormone surges. *J. Clin. Endocrinol. Metab.*1979, 49, 507 – 12.
21. Stricker R et al. Establishment of detailed reference values for luteinizing hormone, follicle stimulating hormone, estradiol, and progesterone during different phases of the menstrual cycle on the Abbott ARCHITECT analyzer". *Clinical Chemistry and Laboratory Medicine.* 2006, 44 (7), 883–7.
22. Simpson ER and McDonald PC. Endocrinology of pregnancy. In: *Textbook of Endocrinology,* Ed.: Williams R.H. 1981, 412 – 22, Saunders Company, Philadelphia.
23. Lipson, SF and Ellison, PT. Comparison of salivary steroid profiles in naturally occurring conception and non-conception cycles. *Hum Reprod,* 1996, 11(10), 2090-96.
24. Choe, JK, Khan-Dawood, FS, and Dawood, MY. Progesterone and estradiol in saliva and plasma during the menstrual cycle. *Am J Obstet Gynecol,* 1983, 147(5), 557-62.
25. Lewis JG. Steroid Analysis in Saliva: An overview. *Clin Biochem Rev* 2006,27 August 139-146.
26. Gameiro, CM, Romão, F, and Castelo-Branco, C. Menopause and aging: Changes in the immune system –A review. *Maturitas,* 2010, 67(4), 316-20.

Please use only the valid version of the Instructions for Use provided with the kit

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